

Cat No: HR1APEA052

For research use only

Overview

Product Name	Akt3 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:500-1,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	60kDakDa
GeneID?Human?	10000
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PKB gamma antibody, PRKBG antibody, Protein Kinase B gamma antibody, RACPK Gamma antibody
Background	RAC-gamma serine/threonine-protein kinase is an enzyme that in humans is encoded by the AKT3 gene. This kinase has been shown to be stimulated by platelet-derived growth factor (PDGF), insulin, and insulin-like growth factor 1 (IGF1).