

## **Bovine Serum Albumin Rabbit pAb** Cat No: HR1APEA054

For research use only

## Overview

Product Name	Bovine Serum Albumin Rabbit pAb
Source	Rabbit
Applications	WB
Species Reactivity	Cow,Human
Recommended Dilutions	WB 1:2,000-5,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	66kDakDa
GeneID?Human?	280717
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BSA antibody, ALB antibody, Albumin antibody, HSA antibody, PRO0883 antibody, Serum Albumin antibody
Background	Bovine serum albumin (also known as BSA or "Fraction V") is a serum albumin protein derived from cows. It is often used as a protein concentration standard in lab experiments. Albumin is a soluble, monomeric protein which comprises about one-half of the blood serum protein.