

# Pan Methyl Lysine Rabbit pAb

Cat No: HR1APEA055

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                   |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Pan Methyl Lysine Rabbit pAb                                                                                                                                                                                                                      |
| Source                | Rabbit                                                                                                                                                                                                                                            |
| Applications          | WB, IHC                                                                                                                                                                                                                                           |
| Species Reactivity    | N/A                                                                                                                                                                                                                                               |
| Recommended Dilutions | WB 1:1,000-2,000 IHC 1:100-200                                                                                                                                                                                                                    |
| Immunogen             |                                                                                                                                                                                                                                                   |
| Species               | Rabbit                                                                                                                                                                                                                                            |
| Storage               | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.                                                                                                                                           |
| Isotype               | IgG                                                                                                                                                                                                                                               |
| Clonality             | Polyclonal                                                                                                                                                                                                                                        |
| Concentration         | 1mg/ml                                                                                                                                                                                                                                            |
| Observed band         | N/AkDa                                                                                                                                                                                                                                            |
| GeneID?Human?         | N/A                                                                                                                                                                                                                                               |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                   |
| Cellular localization |                                                                                                                                                                                                                                                   |
| Alternative Names     | N/A                                                                                                                                                                                                                                               |
| Background            | Many proteins are post translationally modified. Modifications such as phosphorylation, glycosylation, ubiquitination and methylation have been shown to play an important role in the development, physiology and disease of animals and plants. |