

## ZNF499 Rabbit pAb Cat No: HR1APEA057

For research use only

## Overview

Product Name	ZNF499 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:500-1,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	70kDakDa
GeneID?Human?	84878
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Zinc finger protein 499, Zinc finger and BTB domain-containing protein 45
Background	Zinc finger proteins are DNA and RNA binding proteins that contain motifs where amino acids are folded into a single structural unit around a zinc atom. In a standard zinc finger, one zinc atom is bound to two cysteines and two histidines