

## ? Tubulin Rabbit pAb Cat No: HR1APEA058

For research use only

## Overview

Product Name	? Tubulin Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:2,000-5,000 IHC 1:500-1,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	55kDakDa
GeneID?Human?	203068
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Beta 4 tubulin, TBB5, TUBB2, TUBB2A, tubulin beta 2A, beta tubulin, tubulin beta chain
Background	Microtubules are constituent parts of the mitotic apparatus, cilia, flagella, and elements of the cytoskeleton. They consist principally of 2 soluble proteins, alpha- and beta-tubulin, each of about 55,000 Da. Antibodies against beta Tubulin are useful as loading controls for Western Blotting. However it should be noted that levels of ?-Tubulin may not be stable in certain cells. For example, expression of ?-Tubulin in adipose tissue is very low and therefore ?-Tubulin should not be used as loading control for these tissues.