

## Histone H2B(Tri Methyl Lys43) Rabbit pAb Cat No: HR1APEA092

For research use only

## Overview

Product Name	Historia H2P/Tri Mathyl Lyc42\ Pohhit nAh
Product Name	Histone H2B(Tri Methyl Lys43) Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:500-1,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	Antiserum
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	14kDakDa
GeneID?Human?	8349
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	H2B GL105 antibody, H2B.1 antibody, Histone H2B type 1B antibody, HIST1H5BB antibody
Background	Histone H2B is one of the 5 main histone proteins involved in the structure of chromatin in eukaryotic cells.