



Vimentin Rabbit pAb

Cat No: HR1APEA094

For research use only

Overview

Product Name	Vimentin Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:500-1,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	50-57kDa
GeneID?Human?	7431
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FLJ36605 antibody, VIM antibody, VIME antibody
Background	The cytoskeleton consists of three types of cytosolic fibers: microfilaments (actin filaments), intermediate filaments, and microtubules.Vimentin is a type III intermediate filament (IF) protein that is expressed in mesenchymal cells.