

Histone H1(Phospho Thr17) Rabbit pAb Cat No: HR1APEA111

For research use only

Overview

Product Name	Histone H1(Phospho Thr17) Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	Antiserum
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	17-25kDakDa
GeneID?Human?	3005
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	H1 antibody, HIST1H1B antibody, Histone 1 antibody, Histone H1a antibody
Background	Histone H1 is one of the five main histone protein families which are components of chromatin in eukaryotic cells. May play a key role in the control of gene expression during oogenesis and early embryogenesis, presumably through the perturbation of chromatin structure.