

LDHD Rabbit pAb

Cat No: HR1APEA126

For research use only

Overview

Product Name	LDHD Rabbit pAb
Source	Rabbit
Applications	WB
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-5,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	55kDa
GeneID?Human?	197257
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	D Lactate Dehydrogenase antibody, DLD antibody, Lactate dehydrogenase D antibody, MGC57726 antibody
Background	LDHD belongs to the D-isomer specific 2-hydroxyacid dehydrogenase family. D-lactate and D-glycerate dehydrogenase activities have been found in a similar protein in yeast.