

COX IV Rabbit pAb

Cat No: HR1APEA150

For research use only

Overview

Product Name	COX IV Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000 IHC 1:500-1,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	15kDa
GeneID?Human?	1327
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	COX 4, Cox4a, COX4B, COXIV, Cytochrome c oxidase subunit 4 isoform 1 mitochondrial, MGC105470, dj857M17.2, Cytochrome c oxidase subunit IV
Background	The enzyme cytochrome c oxidase or Complex IV, EC 1.9.3.1) is a large transmembrane protein complex found in bacteria and the mitochondrion. It is the last enzyme in the respiratory electron transport chain of mitochondria (or bacteria) located in the mitochondrial (or bacterial) membrane. It receives an electron from each of four cytochrome c molecules, and transfers them to one oxygen molecule, converting molecular oxygen to two molecules of water. In the process, it binds four protons from the inner aqueous phase to make water, and in addition translocates four protons across the membrane, helping to establish a transmembrane difference of proton electrochemical potential that the ATP synthase then uses to synthesize ATP.