

TRIM72/MG53 Rabbit pAb

Cat No: HR1APEA151

For research use only

Overview

Product Name	TRIM72/MG53 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	53kDakDa
GeneID?Human?	493829
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Mitsugumin-53, TRI72, Tripartite motif containing 72, LOC493829
Background	The Tripartite motif(TRIM) family of proteins are characterized by a conserved TRIM domian.Regulates membrane budding and exocytosis. May be involved in the regulation of the mobility of KCNB1-containing endocytic vesicles.