

E-Cadherin Rabbit pAb Cat No: HR1APEA152

For research use only

Overview

Product Name	E-Cadherin Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human
Recommended Dilutions	WB 1:1,000 IHC 1:500-1,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	135kDakDa
GeneID?Human?	999
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CDH-1, Cadherin 1, ECAD, UVO, epithelial calcium dependant adhesion protein, LCAM, CD 234, CD324
Background	E-Cadherin is a calcium-dependent cell-cell adhesion glycoprotein composed of five extracellular cadherin repeats, a transmembrane region, and a highly conserved cytoplasmic tail.