

14-3-3 epsilon Rabbit pAb

Cat No: HR1APEA154

For research use only

Overview

Product Name	14-3-3 epsilon Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Rat,Mouse
Recommended Dilutions	WB 1:1,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	28kDakDa
GeneID?Human?	7531
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	14-3-3 E, 14-3-3 epsilon, FLJ45465, MDCR, YWHAE, protein kinase C inhibitor protein 1, MDS, KCIP 1
Background	The 14-3-3 family of proteins plays a key regulatory role in signal transduction, checkpoint control, apoptotic and nutrient-sensing pathways. 14-3-3 protein epsilon is a protein that in humans is encoded by the YWHAE gene.