



## CNOT4 Polyclonal Antibody

Cat No: HR1AP7516

For research use only

### Overview

Product Name	CNOT4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	64kDa
GenelD?Human?	CNOT4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CNOT4; NOT4; CCR4-NOT transcription complex subunit 4; CCR4-associated factor 4; E3 ubiquitin-protein ligase CNOT4; Potential transcriptional repressor NOT4Hp
Background	CCR4-NOT transcription complex subunit 4(CNOT4) Homo sapiens The protein encoded by this gene is a subunit of the CCR4-NOT complex, a global transcriptional regulator. The encoded protein interacts with CNOT1 and has E3 ubiquitin ligase activity. Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jul 2010],