

AMPK ?1 Rabbit pAb Cat No: HR1APEA160

For research use only

Overview

Product Name	AMPK ?1 Rabbit pAb
Source	Rabbit
Applications	WB
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	63kDakDa
GeneID?Human?	5562
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AAPK1, 5 AMP activated protein kinase alpha 1 catalytic subunit, PRKAA1, SNF1A, MGC33776
Background	AMP-activated protein kinase (AMPK) is highly conserved from yeast to plants and animals and plays a key role in the regulation of energy homeostasis. AMPK is a heterotrimeric complex composed of a catalytic ? subunit and regulatory ? and ? subunits, each of which is encoded by two or three distinct genes (?1, 2; ?1, 2; ?1, 2, 3).