

ERK3 Rabbit pAb

Cat No: HR1APEA163

For research use only

Overview

Product Name	ERK3 Rabbit pAb
Source	Rabbit
Applications	WB
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:500-1,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	83-105kDa
GeneID?Human?	5597
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MAPK6 antibody?p97-MAPK antibody, HsT17250 antibody, MK06 antibody, PRKM6 antibody
Background	Erk3, also known as MAPK6 or p97 MAPK, is almost 50% identical to Erk1/2 at the kinase domain located in its amino-terminal region. Mitogen-activated protein kinase 6 is an enzyme that in humans is encoded by the MAPK6 gene.