

SLUG Rabbit pAb

Cat No: HR1APEA168

For research use only

Overview

Product Name	SLUG Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:500-1,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	30kDa
GeneID?Human?	6591
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MGC10182 antibody, SLUGH1 antibody, SNAI2 antibody, Zinc finger protein SLUG antibody, WS2D antibody
Background	Transcriptional repressor. Involved in the generation and migration of neural crest cells.