



GSK 3? Rabbit pAb

Cat No: HR1APEA169

For research use only

Overview

Product Name	GSK 3? Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Rat, Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	46kDa
GeneID? Human?	2932
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Glycogen synthase kinase 3 beta antibody, GSK3B antibody, Serine/threonine protein kinase GSK3B antibody
Background	Glycogen synthase kinase 3 beta, also known as GSK3B, is an enzyme that in humans is encoded by the GSK3B gene.