

EPG5 Rabbit pAb

Cat No: HR1APEA172

For research use only

Overview

Product Name	EPG5 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	170-292kDa
GeneID?Human?	57724
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Ectopic P granules autophagy protein 5 homolog antibody, HEEW1 antibody, hEPG5 antibody
Background	N/A