

Bestrophin-2 Rabbit pAb Cat No: HR1APEA177

For research use only

Overview

Product Name	Bestrophin-2 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Rat,Mouse
Recommended Dilutions	WB 1:500-1,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	50kDakDa
GeneID?Human?	54831
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BEST 2 antibody, FLJ20132 antibody, Vitelliform macular dystrophy 2 homolog antibody
Background	The Bestrophins are a newly described family of anion channels unrelated in primary sequence to any previously characterized channel proteins. Two different topologies for Bestrophin-1 have been proposed.