

TLR4 Rabbit pAb

Cat No: HR1APEA180

For research use only

Overview

Product Name	TLR4 Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	7099
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ARMD10 antibody, CD284 antibody, Toll like receptor 4 antibody, TOLL antibody
Background	Toll-like receptor 4 is a protein that in humans is encoded by the TLR4 gene.TLR 4 has also been designated as CD284. TLR4 is most abundantly expressed in placenta, and in myelomonocytic subpopulation of the leukocytes.