



TLR1 Rabbit pAb

Cat No: HR1APEA183

For research use only

Overview

Product Name	TLR1 Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	7096
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CD281 antibody, Toll liker receptor 1 antibody, KIAA0012 antibody, LPRS5 antibody
Background	TLR 1 is a member of the Toll-like receptor family (TLR) of pattern recognition receptors of the innate immune system. TLR1 recognizes pathogen-associated molecular pattern with a specificity for gram-positive bacteria. TLR1 has also been designated as CD281.