

TLR8 Rabbit pAb Cat No: HR1APEA186

For research use only

Overview

Product Name	TLR8 Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG UIUEISU
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	51311
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CD288 antibody, MGC119599 antibody, Toll like receptor 8 antibody
Background	Toll-like receptor 8 is a protein that in humans is encoded by the TLR8 gene. TLR8 has also been designated as CD288.