

NMBR (Bombesin Receptor 1) Rabbit pAb

Cat No: HR1APEA187

For research use only

Overview

Product Name	NMBR (Bombesin Receptor 1) Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG DIUEISU
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	50kDakDa
GeneID?Human?	4829
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Neuromedin B receptor antibody, Bombesin Receptor 1 antibody, BB1 antibody
Background	Receptor for neuromedin-B. Belongs to the G-protein coupled receptor 1 family.