

Transferrin Rabbit pAb

Cat No: HR1APEA188

For research use only

Overview

Product Name	Transferrin Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:5,000-10,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	77kDakDa
GeneID?Human?	7018
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Apotransferrin antibody, beta 1 metal binding globulin antibody, serotransferrin antibody, TF antibody
Background	Transferrins are iron-binding blood plasma glycoproteins that control the level of free iron in biological fluids. It is responsible for the transport of iron from sites of absorption and heme degradation to those of storage and utilization.