

CALHM1 Rabbit pAb

Cat No: HR1APEA189

For research use only

Overview

Product Name	CALHM1 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat
Recommended Dilutions	WB 1:500-1,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	40kDaDa
GeneID?Human?	255022
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Calcium Homeostasis modulator 1 antibody, FAM26C antibody, MGC39514 antibody
Background	Calcium homeostasis modulator 1 (CALHM1) is a pore-forming membrane protein that in humans is encoded by the CALHM1 gene.