

## Cav ?3 Rabbit pAb Cat No: HR1APEA191

For research use only

## Overview

Product Name	Cav ?3 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Rat,Mouse
Recommended Dilutions	WB 1:500-1,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	50kDakDa
GeneID?Human?	784
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	3 subunit of voltage-gated Ca2+ channel antibody , Cacnb3 antibody
Background	The beta subunit of voltage-dependent calcium channels contributes to the function of the calcium channel by increasing peak calcium current, shifting the voltage dependencies of activation and inactivation.