

Cav ?4 Rabbit pAb Cat No: HR1APEA195

For research use only

Overview

Product Name	Cav ?4 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	35kDakDa
GeneID?Human?	27092
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	acng4 antibody, TARP gamma-4 antibody, Transmembrane AMPAR regulatory protein gamma-4 antibody
Background	The ? subunit is an integral membrane protein. The ? family consists of at least 8 members, which share a number of common structural features. CaV? subunits inhibit CaV channel activity and modulate its activation and inactivation kinetics.