

Cav ?5 Rabbit pAb

Cat No: HR1APEA196

For research use only

Overview

Product Name	Cav ?5 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:500-1,000 IHC 1:50-100
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	35kDa
GeneID?Human?	10369
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	acng5 antibody, TARP gamma-5 antibody, Transmembrane AMPAR regulatory protein gamma-5 antibody
Background	The ? subunit is an integral membrane protein. The ? family consists of at least 8 members, which share a number of common structural features. CaV? subunits inhibit CaV channel activity and modulate its activation and inactivation kinetics.