

LexA DNA Binding Region Rabbit pAb

Cat No: HR1APEA203

For research use only

Overview

Product Name	LexA DNA Binding Region Rabbit pAb
Source	Rabbit
Applications	WB
Species Reactivity	N/A
Recommended Dilutions	WB 1:10,000-50,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N/A
Background	The LexA protein of Escherichia coli is a transcriptional repressor regulating approximately 20 genes, many of which are involved in DNA repair (SOS response). It binds with variable affinity to single, double, or even triple-operators via its amino-terminal domain that contains three alpha-helices spanning residues.