



## CCDC102A Polyclonal Antibody

Cat No: HR1AP9123

For research use only

### Overview

Product Name	CCDC102A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GenID?Human?	CCDC102A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CCDC102A; Coiled-coil domain-containing protein 102A
Background	CCDC102A (Coiled-Coil Domain Containing 102A) is a Protein Coding gene. An important paralog of this gene is CCDC102B.



**bioelsa**