

PI3 Kinase P85? Rabbit pAb

Cat No: HR1APEA210

For research use only

Overview

Product Name	PI3 Kinase P85? Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:50-100
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	82kDa
GeneID?Human?	5296
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	P85 antibody, p85 beta antibody, PIK3R2 antibody, Phosphoinositide 3 kinase regulatory subunit 2 antibody
Background	PI3 Kinase P85? mediating the association of the p110 catalytic unit to the plasma membrane.