

## c-Raf Rabbit pAb Cat No: HR1APEA215

For research use only

## Overview

Product Name	c-Raf Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	75kDakDa
GeneID?Human?	5894
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	c-Raf proto-oncogene, serine/threonine kinase antibody, EC2.7.11.1 antibody, NS5 antibody, Raf 1 antibody
Background	A-Raf, B-Raf, and c-Raf (Raf-1) are the main effectors recruited by GTP-bound Ras to activate the MEK-MAP kinase pathway. C-Raf is a member of the Raf kinase family of serine/threonine-specific protein kinases, from the TKL (Tyrosine-kinase-like) group of kinases.