

MEK5 Rabbit pAb

Cat No: HR1APEA222

For research use only

Overview

Product Name	MEK5 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	49kDa
GeneID?Human?	5067
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MP2K5 antibody, MAPK kinase kinase 5 antibody, MEK5 antibody, PRKMK 5 antibody
Background	Dual specificity mitogen-activated protein kinase kinase 5 is an enzyme that in humans is encoded by the MAP2K5 gene. The protein encoded by this gene is a dual specificity protein kinase that belongs to the MAP kinase kinase family. This kinase specifically interacts with and activates MAPK7/ERK5.