

HAO1/GOX Rabbit pAb

Cat No: HR1APEA223

For research use only

Overview

Product Name	HAO1/GOX Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Rat, Mouse
Recommended Dilutions	WB 1:2,000-5,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	41kDa
GeneID? Human?	54363
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Glycolate oxidase antibody, GOX1 antibody, HAOX1 antibody, Hydroxyacid oxidase 1 antibody, MGC142225 antibody
Background	Hydroxyacid oxidase (glycolate oxidase) 1 is a protein that in humans is encoded by the HAO1 gene. This gene is one of three related genes that have 2-hydroxyacid oxidase activity yet differ in encoded protein amino acid sequence, tissue expression and substrate preference.