

## Caspase 8 Rabbit pAb

Cat No: HR1APEA225

For research use only

### Overview

Product Name	Caspase 8 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:50-100
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	43,57kDa
GeneID?Human?	841
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ALPS2B antibody, Apoptosis cysteine protease antibody, CAP4 antibody, CED3 antibody, MCH5 antibody
Background	Caspase-8 is a caspase protein, encoded by the CASP8 gene. It most likely acts upon caspase-3. Activated caspase-8 cleaves and activates downstream effector caspases such as caspase-1, -3, -6, and -7.