

STAT1 Rabbit pAb

Cat No: HR1APEA226

For research use only

Overview

Product Name	STAT1 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	84,91kDa
GeneID?Human?	6772
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CANDF7 antibody, ISGF3 antibody, DKFZp686B04100 antibody
Background	STAT1 is a member of the Signal Transducers and Activators of Transcription family of transcription factors. STAT1 is involved in upregulating genes due to a signal by either type I, type II, or type III interferons.