

PPAR ? Rabbit pAb Cat No: HR1APEA228

For research use only

Overview

Product Name	PPAR ? Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:500-1,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	53,57kDakDa
GeneID?Human?	5468
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CIMT1 antibody, NR1C3 antibody, PPARG antibody, PPAR gamma antibody
Background	In the field of molecular biology, the peroxisome proliferator-activated receptors (PPARs) are a group of nuclear receptor proteins that function as transcription factors regulating the expression of genes.