

STAT3 Rabbit pAb

Cat No: HR1APEA234

For research use only

Overview

Product Name	STAT3 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	79,86kDa
GeneID?Human?	6774
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	APRF antibody, HIES antibody, Signal transducer and activator of transcription 3 antibody
Background	Signal transducer and activator of transcription 3, also known as STAT3, is a transcription factor which in humans is encoded by the STAT3 gene. Stat3 is constitutively activated in a number of human tumors and possesses oncogenic potential and anti-apoptotic activities.