

## STAT4 Rabbit pAb Cat No: HR1APEA238

For research use only

## Overview

Product Name	STAT4 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	86kDakDa
GeneID?Human?	6775
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SLEB11 antibody, signal transducer and activator of transcription 4 antibody
Background	STAT4 is a transcription factor belonging to the Signal Transducer and Activator of Transcription protein family. It is required for the development of Th1 cells from naive CD4+ T cells and IFN-? production in response to IL-12.