

TBP Rabbit pAb Cat No: HR1APEA246

For research use only

Overview

Product Name	TBP Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	38kDakDa
GeneID?Human?	6908
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GFT2D antibody, MGC117320 antibody, SCA17 antibody, TATA binding protein antibody, TF2D antibody
Background	The TATA-binding protein (TBP) is a general transcription factor that binds specifically to a DNA sequence called the TATA box. This DNA sequence is found about 30 base pairs upstream of the transcription start site in some eukaryotic gene promoters.