

STAT5b Rabbit pAb

Cat No: HR1APEA247

For research use only

Overview

Product Name	STAT5b Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	90kDa
GeneID?Human?	6777
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	STAT5 antibody, signal transducer and activator of transcription 5b antibody
Background	Signal transducer and activator of transcription 5B is a protein that in humans is encoded by the STAT5B gene. STAT5B orthologs have been identified in most placentals for which complete genome data are available.