

c-Jun Rabbit pAb Cat No: HR1APEA248

For research use only

Overview

| Product Name | c-Jun Rabbit pAb |
|-----------------------|---|
| Source | Rabbit |
| Applications | WB, IHC |
| Species Reactivity | Human, Mouse |
| Recommended Dilutions | WB 1:1,000-2,000 IHC 1:50-100 |
| Immunogen | |
| Species | Rabbit |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Isotype | IgG |
| Clonality | Polyclonal |
| Concentration | 1mg/ml |
| Observed band | 43kDakDa |
| GeneID?Human? | 3725 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | Activator Protein 1 antibody, AP1 antibody, p39 antibody |
| Background | c-Jun is a member of the Jun Family containing c-Jun, JunB and JunD, and is a component of the transcription factor AP-1 (activator protein-1). |