

## c-Jun Rabbit pAb Cat No: HR1APEA248

For research use only

## Overview

Product Name	c-Jun Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human, Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:50-100
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	43kDakDa
GeneID?Human?	3725
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Activator Protein 1 antibody, AP1 antibody, p39 antibody
Background	c-Jun is a member of the Jun Family containing c-Jun, JunB and JunD, and is a component of the transcription factor AP-1 (activator protein-1).