

## S100 Rabbit pAb Cat No: HR1APEA253

For research use only

## Overview

Product Name	S100 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:50-100
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	15kDakDa
GeneID?Human?	6271
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Bpb antibody, NEF antibody, protein S100-A1 antibody
Background	The S100 proteins are a family of low-molecular-weight proteins found in vertebrates and characterized by two calcium-binding sites that have helix-loop-helix ("EF-hand type") conformation.