

Estrogen Receptor ? Rabbit pAb Cat No: HR1APEA256

For research use only

Overview

Product Name	Estrogen Receptor ? Rabbit pAb
Source	Rabbit
Applications	WB
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	66kDakDa
GeneID?Human?	2099
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ER alpha antibody, ESR1 antibody, ESTRR antibody, Major ORF antibody, NR3A1 antibody
Background	Estrogen receptors are a group of proteins found inside cells. Estrogen receptor ? (ER?), a member of the steroid receptor superfamily, contains highly conserved DNA binding and ligand binding domains.