

## Flotillin-2 Rabbit pAb

Cat No: HR1APEA257

For research use only

## Overview

Product Name	Flotillin-2 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:50-100
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	49,58kDakDa
GeneID?Human?	2319
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ECS-1 antibody, ESA antibody, FLOT2 antibody, M17S1 antibody, REG-1 antibody, Reggie-1 antibody
Background	Flotillins belong to a family of lipid raft-associated integral membrane proteins that carry a evolutionarily conserved domain called the prohibitin homology domain. May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.