

# Bcl-Xl Rabbit pAb

Cat No: HR1APEA258

For research use only

## Overview

Product Name	Bcl-Xl Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:50-100
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	30kDa
GeneID?Human?	598
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	B2CL1 antibody, Bcl 2 like 1 antibody, Apoptosis regulator Bcl X antibody
Background	Bcl-xL prevents apoptosis through two different mechanisms: heterodimerization with an apoptotic protein inhibits its apoptotic effect and formation of mitochondrial outer membrane pores help maintain a normal membrane state under stressful conditions.