

Nrf2 Rabbit pAb Cat No: HR1APEA259

For research use only

Overview

Product Name	Nrf2 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	75-100kDakDa
GeneID?Human?	4780
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HEBP1 antibody, Like 2 antibody, NFE2 related factor 2 antibody, Nuclear factor erythroid 2 like 2 antibody
Background	Nuclear factor (erythroid-derived 2)-like 2, also known as NFE2L2 or Nrf2, is a transcription factor that in humans is encoded by the NFE2L2 gene. Small amounts of constitutive nuclear NRF2 maintain cellular homeostasis through regulation of basal expression of antioxidant response genes.