

ASIC1 Rabbit pAb

Cat No: HR1APEA260

For research use only

Overview

Product Name	ASIC1 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:50-100
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	80kDa
GeneID?Human?	41
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ACCN2 antibody, BnaC2 antibody, Brain sodium channel 2 antibody
Background	Amiloride-sensitive cation channel 2, neuronal is a protein that in humans is encoded by the ACCN2 gene.